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Analysis of Small Amounts of Harmful Air Pollutants

Fukuo Sagawa・Masayuki Taniguchi・Shigeo Tachibana・Haruo Yoshida

Recent analytical techniques used for monitoring air pollution by harmful substances were described. The substances covered were 11 volatile organic compounds (VOCs) and two metallic elements (Ni, As) which were selected as pollutants requiring urgent control in Japan. The analytical techniques used were mainly the canister method for sampling, the GC/MS method for VOCs analysis and the ICP method for elemental analysis. Satisfactory analytical capabilities at Kobe Steel Group were also mentioned.

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Table 2 Environmental standards of VOCs in atmosphere

VOCs	Standard Value*)	Analytical Method
Benzene	0.003 mg/m ³	Sampling : canister or collecting tube method Analysis : GC/MS method or equivalents
Trichloroethylene	0.2 mg/m ³	
Tetrachloroethylene	0.2 mg/m ³	

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
Fig. 1 Sampling methods by use of canisters

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Table 3 Analytical methods for harmful atmospheric pollutants

Pollutant	Sampling Method	Analytical Method
Benzene Trichloroethane Tetrachloroethane Dichloromethane	Collection of air sample by evacuated stainless steel canister	Concentration at low temperature followed by GC/MS analysis
Vinyl Chloride Acrylonitrile Chloroform 1,2-dichloromethane 1,3-butadiene	Suction of air sample by a pump through collecting tube packed with suitable adsorbents	Desorption by heating followed by GC/MS analysis
Benzene Trichloroethane Tetrachloroethane Dichloromethane Acrylonitrile	Suction of air sample by a pump through collecting tube packed with suitable adsorbents	Solvent extraction followed by GC/MS analysis
Acetaldehyde Formaldehyde	Suction of air sample by a pump through collection tube packed with a chemical reactant	Solvent extraction followed by GC/MS, HPLC or GC (FTD) analysis
Nickel & Its Compounds Arsenic & Its Compounds	Collection of particulate substances on a filter by use of high-volume or low-volume air sample	Acid extraction or hydride generation followed by AAS, ICP/AES or ICP/MS analysis

HPLC : High Performance Liquid Chromatography
AAS : Atomic Absorption Spectrometry
AES : Atomic Emission Spectroscopy

FTD : Flame Thermionic Detector
ICP : Inductively Coupled Plasma

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Method	Atomic Absorption (Flame)	Atomic Absorption (Flameless)	ICP Atomic Emission (ICP-AES)	ICP Mass Spectrometry (ICP-MS)
Determination Range ($\mu\text{g}/\text{ml}$) ^{*1}	0.3 ~ 6	0.005 ~ 0.05	0.04 ~ 2	0.001 ~ 0.1
Preparation of Sample Solution ^{*2}	A, B, C, D, F	A, B, C, D, F	A, B, C, D, F	B ^{*3}
Possible Elements for Simultaneous Determination	Cu, Zn, Pb, Cd, Mn, Fe, Co, etc.	Cu, Zn, Pb, Cd, Mn, Fe, Co, Be, etc.	Cu, Zn, Pb, Cd, Mn, Fe, Co, Be, etc.	Cu, Zn, Pb, Cd, Mn, Fe, etc.

*1 An example of range shown as concentration in sample solution

*2 A : Dissolution by HF·HNO₃·HC/O₄ solution
B : Dissolution in pressurized vessel
C : Dissolution by HF·H₂O₂ solution
D : Dissolution by HNO₃·HC/ solution
F : Solvent extraction

*3 Other methods can be applicable but B method is the basic one.

第 5 表 ヒ素分析の各種方法
Table 5 Methods for arsenic analysis

Method	Hydride Generation & Atomic Absorption	Hydride Generation & ICP-AES	ICP-MS (Direct Analysis)
Determination Range ($\mu\text{g}/\text{ml}$) ^{*1}	0.005 ~ 0.05	0.005 ~ 0.05	0.001 ~ 0.1
Preparation of Sample Solution ^{*2}	B, E	B, E	B
Possible Elements for Simultaneous Determination	Se	Se	Cu, Zn, Pb, Cd, Mn, Co, etc.

*1 An example of range shown as concentration in sample solution

*2 B : Dissolution in pressurized vessel
E : Dissolution by HNO₃·H₂SO₄ solution

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Table 6 Apparatus used for analysis of small amounts of harmful air pollutants (in Kobe Steel Group)

Substance Analyzed	Main Apparatus Used
VOCs	GC/MS (MS: Quadrupole), HPLC, GC (Detector: FTD)
Ni, Cr and Their Compounds	AAS (Flame & Flameless Methods), ICP/AES, ICP/MS (MS: Quadrupole)
Asbestos	EDX/SEM, TEM
Dioxin	GC/MS (MS: Magnetic Sector)

EDX: Energy Dispersive X-ray Spectrometer
TEM: Transmission Electron Microscope

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
Fig. 4 An analytical result (chromatogram) of mixed standard VOCs obtained by GC/MS method

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