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Thin Film Analysis using a High Resolution RBS System

笹川 薫  
Kaoru Sasakawa

High resolution RBS systems have been applied for the analysis of ultra-thin films in various industrial fields. In this paper, recent ultra-thin film application examples are given. These include high permittivity films for next-generation gate dielectric films, distortion analysis or damage in crystal surfaces or interfaces, and hydrogen analysis.

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Table 1 Analytical needs of ultra thin films by high-resolution RBS

Product	Material	Notable layer, Thickness	Evaluation item	Related film character
Silicon device	High-k gate dielectric	Interface layer, 1-5 nm	Depth profile	Equivalent oxide thickness
	SiON gate dielectric	Nitride layer, 1-3 nm	Depth profile of nitrogen	Reliability
	SiO <sub>2</sub> gate dielectric	Within SiO <sub>2</sub> and interface	Interface strain, areal density	Reliability
	Si/SiGe	Epitaxial Si	Depth profile of strain, composition of SiGe	Mobility
	Shallow junction	Doped layer, 10-30 nm	Ion dose, damage, position of implanted atom	Resistance
	SiONO	Each interface, 20-30 nm	Depth profile, areal density of each layer	Capacitance
	Barrier TiN, TaN	Within film, 3-10 nm	Depth profile, areal density	Reliability
	Wafer	Near surface	Depth distribution of damage and defect	Reliability
Hard disk drive	Head (GMR,TMR)	Each layer, 1-3 nm	Depth profile and areal density	Magneto resistance
	Magnetic film	Each layer, 1-50 nm	Depth profile and areal density of each layer	Recording property
	Overcoat	Within film, 1-5 nm	Depth distribution of hydrogen, areal density	Tribology, wear resistance
	Lubricant	Within layer, 1-3 nm	Areal density, coverage	Tribology
Optical device	LED,LD	Active layer (MQW <sup>*1</sup> ), 1-3 nm	Composition and areal density	Emission property
	Optical crystal	Near surface	Depth distribution of damage and defect	Optical transmittance

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